

# LaserLED 16 Infrared Light Torch

## Directions for use:

The light can be used on any part of the body for any length of time with complete safety except for the eyes or face where we recommend that you use protective eye-wear suitable for the invisible Infrared light frequencies.

1. **To treat for pain or any other condition** simply place the Light directly on to the skin and hold in position for anywhere between 30 seconds to 2 minutes.

It is possible that a longer duration of treatment time may be of additional benefit. However, research has indicated that the millions of cells being treated in the immediate vicinity of the application may not absorb any more benefit with sustained use of a single continuous application.

2. **Repeat treatments** after a 2 hour interval may be applied if necessary.

## THE MERIDIAN SYSTEM:

Acupuncturists discovered that single frequency light could activate acupressure points. Pulse light could stimulate it; continuous light could sedate the acupuncture points. But they also discovered that light applied to a *meridian* end-point can actually be traced flowing through the meridian to the organ acupuncture points. The meridian system is a useful pathway for getting light deeper into the body, so if you are treating things like asthma in the lungs, there is an alternate method of getting light into the lungs.

Excellent, simple to follow  
Acupuncture Charts are included with the...

## **'Search N Stim' Electronic Acupuncture Point DETECTOR**

which can be purchased as an  
Optional extra from our website for  
only \$125 incl. GST & p&p within  
Australia.

<http://tinyurl.com/ecdsf>

Some Important Acupuncture Charts &  
Conditions used in conjunction with the

## LaserLED 16 Infrared Light Torch:

- Pain in the arms or elbows
- Pain in the hands
- Pain in the legs
- Pain in the shoulder
- Pain in the neck
- Toothache
- Headache, migraine
- Menstrual pain
- Nausea – motion sickness
- Insomnia
- Sexual vitality
- Stress
- Smoking Cessation
- Appetite suppression for weight control

**Download illustrated Acupuncture Points -**  
<http://www.acuxo.com/index.asp>

## PENETRATING THROUGH THE

**BLOOD STREAM:** You can even get light into the blood stream. One of the best ways is through your belly button, because the aorta artery is behind the navel. So if you insert the Light there for 20 minutes, every drop of blood in the body will pass in front of the light, increasing the activity of your white cells, red cells, B-cells and T-cells, so you can boost your whole immune system.

**Another method** to achieve the same result is to place the Array over the Radial & Ulna arteries for 3-4 minutes each wrist. (The centre of the wrist flex creases.) This is particularly valuable for **Body Builders & Weight Lifters** to utilise the immediate effect of Nitric Oxide production!

**No claims** are made that the LaserLED Infrared Light will treat any medical condition.

**Please consult** your medical adviser before using the Array if you have any questions relating to your particular medical needs

**Use care in the treatment of pregnant women.** (Unless you are a qualified Acupuncturist do not apply this Light to points that coincide with acupuncture points affecting pregnancy particularly in the lower abdominal area.)

**To date there have been no reports whatsoever of any negative effects using an LED Infrared Light device.**

## Good Vibrations

PO Box 4174 Eight Mile Plains  
Brisbane Queensland 4113  
Australia  
Ph: 617 3423 0066  
Free Call: 1800 307 288

[www.goodvibrations-energiser.com](http://www.goodvibrations-energiser.com)  
[ledinfo@goodvibrations-energiser.com](mailto:ledinfo@goodvibrations-energiser.com)

April 2006